

Test changes

Test name: Clostridium difficile toxin PCR

Test #: 8288/LAB8288

Notification date: December 20, 2019

Effective date: January 14, 2020

EXPLANATION

On January 14, 2020, Allina Health Laboratory will implement a multi-step algorithm for *C. difficile* testing

DETAILS

Effective January 14, 2020 Allina Health Laboratory will be instituting a multi-step algorithm for testing and reporting *C. difficile* results. Positive *C. difficile* toxin PCR tests will reflex additional testing, at an additional charge, to confirm that the patient has existing *C. difficile* disease as opposed to being colonized or a carrier.


The PCR test detects the gene encoding for *C. difficile* toxin but does not indicate whether the toxins are actively being produced. The reflexed EIA (enzyme immunoassay) test will detect GDH antigen, indicating the presence of the organism. In addition, it will detect the presence of the Toxin A/B.

By using a testing algorithm that utilizes PCR as the first layer of testing the lab is able to screen out most of the negative samples. If a sample is PCR positive it is reflexed to a less sensitive EIA test to screen out possible carriers, a positive EIA with antigen and toxin detected are considered a true *C. difficile* infection.

Note: Questions regarding interpretation of this algorithm should be directed to Infection Prevention or an Infectious Disease practitioner. The lab cannot interpret patient specific result for this algorithm.

The testing results will be accompanied by an interpretation and comment when indicated:

Positive for toxigenic *C. difficile*

ID STOOL	
CTA INTERPRETATION	Positive for to... !
C. DIFFICILE GDH A...	Positive
C. DIFFICILE TOXIN...	Positive
CLOSTRIDIUM DIFFIC...	

Negative for toxigenic C. difficile

12/12/2019 1427	
ID STOOL	
CTA INTERPRETATION	Negative for to...
C. DIFFICILE GDH A...	Negative
C. DIFFICILE TOXIN...	Negative
CLOSTRIDIUM DIFFIC...	

Suggests no Active C. difficile disease

Comment: C. difficile present but toxin was not detected. This could be a non-toxigenic strain, a carrier state, or toxin level too low for detection. The result should be used in conjunction with clinical history when deciding on additional testing and/or possible therapy. Consider ID consultation.

ID STOOL	
CTA INTERPRETATION	Suggests no ac... *
C. DIFFICILE GDH A...	Positive
C. DIFFICILE TOXIN...	Negative
CLOSTRIDIUM DIFFIC...	

de-

Invalid

Comment: Invalid result. Consider recollection if clinical suspicion of infection still exists.

ID STOOL	
CTA INTERPRETATION	Invalid *
C. DIFFICILE GDH A...	
C. DIFFICILE TOXIN...	
CLOSTRIDIUM DIFFIC...	

Printed report

CLOSTRIDIUM DIFFICILE TOXIN P (Final result)			
ID:	19CL-347M0001	Collected:	12/13/2019 1117
Source:	Stool	Verified On:	12/13/2019 1216
Resulting Lab:	ANWLCR		
Component		Value	Ref. Interval
CLOSTRIDIUM DIFFICILE PCR		Positive	
Positive C. difficile PCR results reflexed to EIA testing for confirmation.			
PRESUMPTIVE NAP1 STRAIN		Negative	
Comments: The NAP1 (027 or BI) strain is a hypervirulent strain. Detection may be useful for epidemiological purposes.			
Toxin Activity (Final result)			
ID:	19CL-347M0001	Collected:	12/13/2019 1117
Source:	Stool	Verified On:	12/13/2019 1216
Resulting Lab:	ANWLCR		
Component		Value	Ref. Interval
CTA INTERPRETATION		Suggests no active C. difficile disease (A)	Negative for toxigenic C. difficile
C. difficile present but toxin was not detected. This could be a non-toxigenic strain, a carrier state, or toxin level too low for detection. The result should be used in conjunction with clinical history when deciding on additional testing and/or possible therapy. Consider ID consultation.			
C. DIFFICILE GDH ANTIGEN EIA		Positive	
C. DIFFICILE TOXIN A & B EIA		Negative	

CPT Code for EIA Testing: 87324 and 87449