

**FIBRINOGEN, QUANTITATIVE**  
101/85384.0

**Quick Reference**

Due to implementation of a new reagent lot number, effective June 21, 2022, the Fibrinogen, quantitative reference range will be changing.

Changing to a new lot of coagulation assay reagent requires thorough validation testing to establish a new mean normal reference value. This process occurs approximately every 12 months when a new lot of reagent is made available by the instrumentation reagent vendor.

**Details**

<b>Current (mg/dL)</b>	<b>NEW (mg/dL)</b>
233 - 384	240 - 410

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department