

FIBRINOGEN, QUANTITATIVE
101/85384.0

QUICK REFERENCE

Changing to a new lot of coagulation assay reagent requires thorough validation testing to establish a new mean normal reference value. This process occurs approximately every 12 months when a new lot of reagent has been made available by the instrumentation reagent vendor.

DETAILS

Effective May 18, 2021 the reference range for the Fibrinogen, quantitative assay (101/85384.0) will be changing.

Current	NEW
200 - 450 mg/dL	233 - 384 mg/dL

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department