

PROTEIN ELP, RANDOM URINE
3218/LAB3218

QUICK REFERENCE

On February 21, 2022, there will be a change to the specimen requirements for the Protein ELP, random urine (3218/LAB3218) assay.

DETAILS

	Current	NEW
Collection container	Screw cup plastic container (non-sterile)	Screw cap plastic container (sterile)
Volume	8mL	20mL
Minimum volume	5mL	10mL
Processing instructions	1. Transfer an 8mL aliquot to a Urine aliquot transport vial/tube – 10mL (AHL) 2. Label the tube as urine	<i>None</i>
Transport container	Urine aliquot transport vial/tube – 10mL (AHL)	Screw cap plastic container (sterile)

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department