

Test update

Test name: Haptoglobin
Test #: 127
Notification date: May 28, 2020
Effective date: June 2, 2020

QUICK REFERENCE

Effective June 2, 2020, the methodology for the Haptoglobin (127) assay will change from Rate nephelometry to Immunoturbidimetric.

DETAILS

| | Current | NEW |
|---------------------------------------|-------------------|------------------------|
| Specimen type (preferred) | Serum | Lithium heparin plasma |
| Collection container | Gold SST | Lt green PST |
| Alternate specimen type | EDTA plasma | Serum |
| Alternate collection container | Red, lavender | Gold SST |
| Performing lab | AHL – Immunology | AHL – Chemistry |
| Methodology | Rate Nephelometry | Immunoturbidimetric |
| Reference range | 30 – 200mg/dL | 30 – 215 mg/dL |