

**ESTRONE, SERUM**  
13327/LAB13327

**QUICK REFERENCE**

LabCorp has shared that on April 4, 2022, multiple changes to the Estrone, serum assay (13327/LAB13327) will go into effect.

**DETAILS**

**SPECIMEN TYPE/COLLECTION CONTAINER**

Current: Serum/Gold or red  
EDTA plasma/Lavender

NEW: Serum/Gold *only*

**SPECIMEN VOLUME**

Current: 1.0 mL (minimum 0.6mL)

NEW: 0.6 mL (minimum 0.4 mL).

*Submission of the minimum volume does not allow for repeat testing*

**TRANSPORT/STABILITY**

Current: Refrigerated

NEW: Ambient (*preferred*) – 14 days

Refrigerated – 14 days

Frozen – 14 days

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department

**REFERENCE RANGE**

**CURRENT**

	Male (pg/mL)	Female (pg/mL)
Tanner Stage I	5–17	4–29
Tanner Stage II	10–25	10–33
Tanner Stage III	15–25	15–43
Tanner Stage IV	15–45	16–77
Tanner Stage V	20–45	29–105
Adults	15–65	
Follicular phase:		39–132
Periovulatory:		58–256
Luteal phase:		54–179
Pregnancy		
1st trimester:		247–2774
2nd trimester:		569–5781
Postmenopausal		
with ERT:		51–488
without ERT:		31–100

**NEW**

Age (y)	Male (pg/mL)	Female (pg/mL)
Tanner Stage I	5–17	4–29
Tanner Stage II	10–25	10–33
Tanner Stage III	15–25	15–43
Tanner Stage IV	15–45	16–77
Tanner Stage V	20–45	29–105
Adult	0 - 174	
Premenopausal		27 - 231
Menstrual cycle (1 – 10 days)		19 - 149
Menstrual cycle (11 – 20 days)		32 - 176
Menstrual cycle (21 – 30 days)		37 - 200
Postmenopausal		0 - 125

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department