

New test

Test name:	Coccidioides Ab, reflexive panel
Test #:	994/LAB994
Notification date:	November 17, 2020
Effective date:	Immediately

Effective immediately, a new Coccidioides Ab, reflexive panel, referred to LabCorp, is available. This test is ordered as a Miscellaneous Sendout (994/LAB994).

Test name:	Coccidioides Ab, reflexive panel
Test number:	994
Excellian order number:	LAB994

Abbreviation:	MSO
Tests included:	<ul style="list-style-type: none">• Coccidioides antibody IgG and IgM by ELISA are used to screen for Coccidioides antibodies.• If the ELISA testing is equivocal or positive for IgG and/or IgM, then <i>Coccidioides immitis</i> antibodies by Immunodiffusion, and Coccidioides titer will be added.

Useful for: Preferred serology test to detect Coccidioidomycosis (Valley fever)

Collection container: [Gold serum separator \(SST\) tube](#)



Volume:	2.0 mL
Minimum volume:	0.6 mL
Processing instructions:	Spin
Transport container:	Gold serum separator (SST) tube
Transport and stability:	Refrigerated (<i>preferred</i>) – 2 weeks Ambient – 48 hours Frozen – 1 year