

Allergen Listing

ImmunoCAP IgE

Test name	test #	Test name	test #	Test name	test #
Adult food allergy profile	12452	Cheese, cheddar type	13453	Fenugreek seed, IgE	994
Almond (F20) IgE	4223	Cheese, mold type, IgE	994	Ferret epithelium, IgE	994
Alpha lactalbumin	13452	Chestnut, sweet (F299), IgE	12444	Flounder, IgE	994
Alpha-gal meat allergy panel	994	Chick pea, IgE	13464	Fusarium moniliforme, IgE	12445
Alternaria alternata (M6), IgE	4195	Chicken feathers, IgE	994	Garlic, IgE	13439
Anchovy, IgE	994	Chicken meat (F83), IgE	4262	Gelatin, IgE	994
Aniseed, IgE	994	Childhood allergy march profile	12451	Gelatin, porcine, IgE	994
Annatto seed, IgE	994	Chili Pepper, IgE	994	Gerbil epithelium, IgE	994
Apple (F49), IgE	4898	Chocolate (F93), IgE	4236	Giant ragweed (W3), IgE	4242
Apricot, IgE	994	Cinnamon, IgE	994	Ginger, IgE	994
Aspergillus fumigatus IgE	4206	Cladosporium herbarum (M2), IgE	4198	Gluten (F79), IgE	4255
Aureobasidi pullulans	13469	Clam (F207), IgE	4228	Goat epithelium, IgE	994
Avocado, IgE	13456	Cloves, IgE	994	Goldenrod, IgE	994
Banana (F92), IgE	4904	Cockatiel droppings, IgE	994	Goose feathers, IgE	994
Barley IgE	4906	Cockatiel feathers, IgE	994	Grape, IgE	994
Bass, Black, IgE	994	Cocklebur, IgE	13473	Grapefruit, IgE	994
Beef (F27), IgE	8562	Cockroach (I6), IgE	4249	Green bell pepper, IgE	994
Bermuda grass, IgE	13434	Coconut, IgE	13438	Green pea (F12), IgE	8572
Beta-lactoglobulin, IgE	994	Codfish (F3), IgE	4221	Green string bean (F315), IgE	8568
Bing Cherry, IgE	994	Coffee, IgE	994	Grouper (Serranidae family), IgE	994
Birch (T3), IgE	4213	Common ragweed (W1), IgE	4260	Guava, IgE	994
Black bean, IgE	994	Coriander/Cilantro, IgE	994	Guinea Pig epithelium, IgE	994
Black Pepper, IgE	994	Corn (F8), IgE	4229	Gum Guar, IgE	994
Blackberry, IgE	994	Cottonseed, IgE	994	H halodes (M8), IgE	4231
Blueberry, IgE	13462	Cottonwood (T14), IgE	4244	Haddock, IgE	994
Brazil nut (F18) IgE	4241	Cow dander, IgE	994	Halibut, IgE	13463
Brazil nut IgE with component	994	Crab (F23), IgE	4212	Hamster epithelium, IgE	994
Broccoli, IgE	994	Cranberry, IgE	994	Hazelnut (F17), IgE	4230
Buckwheat, IgE	994	Crayfish, Freshwater, IgE	994	Hazelnut, IgE with component reflex	994
Budgerigar droppings, IgE	994	Cucumber, IgE	994	Hickory, white, IgE	13471
Budgerigar feather, IgE	994	Cumin, IgE	994	Honey bee (I1), IgE	4243
Candida albicans (M5), IgE	4248	Date, IgE	994	Honey, IgE	994
Carageenan gum, IgE	994	Dermatophagoides farinae, IgE	4194	Hop (food), IgE	994
Cardamom seed, IgE	994	Dermatophagoides pteronyssinus, IgE	4207	Hornet, white face, IgE	13465
Carmine red dye, IgE	994	Dill, IgE	994	Hornet, yellow, IgE	13468
Carrot, IgE	13440	Dockweed, Yellow, IgE	994	Horse dander, IgE	13447
Casein (F78), IgE	4239	Dog dander (E5), IgE	4196	House dust (Greer), IgE	4219
Cashew nut (F202), IgE	4224	Duck feathers, IgE	994	House dust HS (H2), IgE	4214
Cashew nut, IgE with component reflex	994	Egg mix, yolk & white (F245), IgE	4240	Immunoglobulin E (IgE)	137
Cat dander (E1), IgE	4193	Egg white (F1), IgE	4205	Insulin, human, IgE	994
Catfish, IgE	994	Egg yolk (F75), IgE	4215	Ispaghula, IgE	994
Cauliflower, IgE	994	Eggplant, IgE	994	Jalapeno pepper, IgE	994
Cedar, IgE	994	Elm (T8), IgE	4218	Johnson grass (G10), IgE	4256
Cedar, Mountain, IgE	994	English plantain (W9), IgE	4251	June (Kentucky) grass, IgE	4261
Celery, IgE	994	Epicoccum purpur, IgE	13445	Kidney bean, IgE	13461

Allergen Listing

ImmunoCAP IgE

Test name	test #	Test name	test #	Test name	test #
Kiwi fruit, IgE	13454	Pear, IgE	994	Squid, IgE	994
Kochia, IgE	994	Pecan nut (F201), IgE	4225	Stachybotrys atra, IgE	994
Lamb's quarter (W10), IgE	4200	Penicillin G, IgE	994	Strawberry (F44), IgE	4237
Latex (K82), IgE	4202	Penicillin V, IgE	994	Sunflower seed, IgE	12449
Lemon, IgE	994	Penicillium chrysogenum (M1), IgE	4208	Sweet Chestnut, IgE	994
Lentil, IgE	13460	Perch, Ocean, IgE	994	Sweet potato, IgE	994
Lima Bean, IgE	994	Perennial rye grass, IgE	4222	Swine epithelium, IgE	994
Lime, IgE	994	Phoma betae, IgE	13444	Swordfish, IgE	994
Linseed, IgE	994	Pigweed, common, IgE	13436	Sycamore, Maple leaf, IgE	13472
Lobster (F80), IgE	4220	Pine nut, pignoles, IgE	12447	Tapioca, IgE	994
Macadamia nut, IgE	13446	Pine, White, IgE	994	Tea, IgE	994
Mackerel, IgE	994	Pineapple, IgE	13457	Tilapia, IgE	13443
Mango, IgE	994	Pistachio nut (F203), IgE	7054	Timothy grass (G6), IgE	4199
Maple (Box elder) (T1), IgE	4216	Plum, IgE	994	Tomato (F25), IgE	4238
Meadow fescue, IgE	994	Pollock, White, IgE	994	Tree nut profile, w/ component reflex	994
Melon, IgE	13455	Pomegranate, IgE	994	Trout, IgE	994
Mesquite, IgE	994	Poppy seed, IgE	994	Tuna (F40), IgE	4227
Milk (F2), IgE	4201	Pork (F26), IgE	8573	Turkey, IgE	994
Millet, IgE	994	Potato, white, IgE	13437	Turmeric, IgE	994
Mold panel	994	Pumpkin Seed, IgE	994	Vanilla, IgE	994
Mouse epithelium, IgE	13475	Pumpkin, IgE	994	Walleye pike, IgE	13451
Mouse serum, IgE	13449	Rabbit epithelium, IgE	994	Walnut (F256), IgE	4217
Mouse urine, IgE	13448	Ragweed, western, IgE	994	Walnut tree (T10), IgE	6795
Mucor racemosus (M4), IgE	4258	Raspberry, IgE	994	Walnut, IgE with component reflex	994
Mugwort, IgE	12446	Red beet, IgE	994	Wasp, paper, IgE	13467
Mulberry, White, IgE	994	Red snapper, IgE	994	Watermelon, IgE	994
Mushroom, IgE	994	Red top (Bent) grass, IgE	4264	Wheat (F4), IgE	4209
Mussel, IgE	13450	Respiratory disease allergen profile	12450	White ash (T15), IgE	6725
Mustard, IgE	994	Rhizopus nigricans, IgE	12448	White bean, IgE	13442
Nettle, IgE	13474	Rice (F9), IgE	4235	White oak (T7), IgE	4197
Nut (food) panel # 2	994	Rough marsh elder (W16), IgE	6872	Willow, IgE	13470
Nutmeg, IgE	994	Russian thistle (W11), IgE	6880	Yeast, IgE	994
Oat (F7), IgE	4263	Rye, IgE	13441	Yellow Jacket, IgE	13466
Octopus, IgE	994	Safflower, IgE	994		
Onion, IgE	994	Salmon (F41), IgE	4226		
Orange, IgE	13435	Sardine, IgE	994		
Orchard grass (G3), IgE	6383	Scallop (F338), IgE	8139		
Ovalbumin, IgE	13458	Seasonal inhalants allergen profile	994		
Ovomucoid, IgE	13459	Sesame seed (F10), IgE	8327		
Oyster (F290), IgE	4252	Setomelanomma rostrata, IgE	994		
Papaya food, IgE	994	Sheep epithelium, IgE	994		
Paprika, IgE	994	Sheep sorrel, IgE	13642		
Parsley, IgE	994	Shrimp (F24), IgE	4211		
Peach, IgE	994	Soybean (F14), IgE	4210		
Peanut (F13), IgE	4204	Spinach, IgE	994		
Peanut components allergen profile	12468	Spruce, IgE	994		

Tests in bold are performed at the AHL Central Laboratory

July 2019