

Test name: Platelet autoantibody profile, whole blood

Test #: 994/LAB994

Notification date: December 24, 2019

Effective date: December 24, 2019

The following new test, [Platelet autoantibody profile, whole blood](#), is available for order effective immediately.

Test name:	Platelet autoantibody profile, whole blood
Test number:	994
Excellian order:	LAB994
Abbreviation	MSO
Useful for:	Detect the presence of an autoantibody that has been associated with autoimmune thrombocytopenic purpura (AITP)
Test includes	Platelet-associated antibodies to platelet-specific glycoproteins IIb/IIIa, Ib/IX, and Ia/IIa
Collect:	12 mL EDTA whole blood
Container:	Lavender (EDTA)
Transport/Stability:	Maintain specimen at room temperature – 72 hours
Alternate Names:	Direct Platelet Antibodies Platelet-associated Immunoglobulin
Performing Lab:	LabCorp (014102): R-NX
Expected TAT:	2 - 5 days
Reference Range:	Negative
Method:	Solid-phase enzyme-linked immunoassay after elution of antibodies from platelets
CPT Code(s):	86022

Interfaced clients:

Please contact your account representative to arrange for build and testing of this assay if this is a test that is utilized at your practice.