

TROPONIN T
TEST ID: 13395/LAB13395

QUICK REFERENCE

LabCorp has announced changes to the Troponin T (13395/LAB13395) assay that will go into effect on February 1, 2021.

DETAILS

	Current	<i>NEW</i>
Collection container	Gold serum separator (SST) tube	Lt green plasma separator (PST)
Alternate collection container	None	Dk green Lithium heparin (Li hep), no gel
Causes for rejection	None	Sample other than heparin plasma
Reference range	<0.011 ng/mL	Male: 0-22 ng/L Female: 0-14 ng/L

QUESTIONS: Contact your Allina Health Laboratory account representative, or our Client Services department